



BreatheSafe® Bite Block Technical Specifications

Material	Medical Grade TPE (Thermoplastic Elastomer) - Latex Free	
Size 1 (BB01-01)	Aqua green, white handle	Small jaw and dentition (7 yrs and over)
Size 2 (BB01-02)	Aqua green, marine green handle	Medium jaw and dentition
Size 3 (BB01-03)	Aqua green, orange handle	Large jaw and dentition
Directions for use:	For Single-Use application only	
Packaging:	Sterile Supplied in individual sealed pouches	Transit Box 200 pcs Carton 50 pcs
Force Resistance:	Over 1500N – exceeds maximum human bite force (1240N)	
Marking:	CE 0120	
Handling & Storage:	Normal hospital protocol for disposal of clinical waste should be followed for disposal of this product. Store between -40°C and +60°C, in dry conditions and away from direct sunlight.	



STERILE R



Read accompanying instructions

BreatheSafe® Bite Block

BreatheSafe® is a registered trademark of Owen Greenings & Mumford Ltd

OGM Ltd
Units 1-4 Oxford Industrial Park
Mead Road
Yarnton
OXFORD OX5 1QU

T +44 (0)1865 844300
E breathesafe@ogm.uk.com
W breathesafe.uk.com



Airway care from start to finish

